

Sponsors

We thank the following sponsors:

Platinum

Bayer Animal Health
National Pork Board
Pfizer Animal Health

Silver

Boehringer Ingelheim Vetmedica, Inc.

Bronze

Cargill
Merck Animal Health
Novartis Animal Health

Copper

AgStar Financial Services
Elanco Animal Health
IDEXX
Newport Laboratories
PIC USA
PRRS CAP

University of Minnesota Institutional Partners

College of Veterinary Medicine
University of Minnesota Extension
College of Food, Agriculture and Natural Resources Sciences

Formatting

Tina Smith Graphics
www.tinasmithgraphics.com

CD-ROM

David Brown
www.davidbrown.us

Logo Design

Ruth Cronje, and Jan Swanson;
based on the original design by Dr. Robert Dunlop

The University of Minnesota is committed to the policy that all persons shall have equal access to its programs, facilities, and employment without regard to race, color, creed, religion, national origin, sex, age, marital status, disability, public assistance status, or sexual orientation.

Detecting Highly Virulent Chinese-Type Porcine Reproductive and Respiratory Syndrome Virus (H-PRRSV) by a Field PCR Platform – POCKIT

Te Yu Chung, Hsiao Fen Grace Chang
GeneReach Biotechnology Corporation, Taiwan, R.O.C.

A highly sensitive field diagnostic platform is an essential tool for farmers. It allows quick decision making to avoid significant economic impact. POCKIT is a qualitative PCR detection system based on insulated isothermal PCR (iiPCR). The PCR reaction is driven by the Rayleigh-Bénard convection principle. A natural thermal convection phenomenon is induced inside a capillary tube with a single heating source applied at the bottom, which results in a temperature gradient between 60°C to 95°C inside the tube. The convection will repeatedly circulate the reagents through different temperature zones for the three steps of PCR –denaturing, annealing, and extension respectively.

A field test was conducted to detect the highly virulent Chinese-type PRRSV (H-PRRSV) by POCKIT. A set of specific primers and a FAM labeled TaqMan probe were designed to detect the H-PRRSV but not the normal strains including Asian origin, NA origin, and EU origin. The RNA was prepared by a 15-minute spin-column protocol. The result showed that the POCKIT system was capable of detecting the field samples efficiently. When confirmed and quantified by the real-time PCR, the titers of the field samples were in the range between 10^3 to 10^5 viral particles/ μ l of extracted samples. A normal PRRSV BSL vaccine strain was introduced as a specificity control and the result was negative accordingly.

The POCKIT system is specifically designed for field diagnosis though it can be used for reference labs for screening tests. It comes in as

a portable hard-case carry-on suitcase with a mini-centrifuge and two micro pipettes. The total run time from sample to result is less than 1.5 hours. With its high sensitivity and specificity, easy-of-use and short-turn-around-time, POCKIT can provide farmers and field consultants a powerful tool for veterinary disease management.